

YTA SYBR Green qPCR MasterMix 2X

Cat No :YT2551

Size: 1 ml

For Research Use Only

Description: YTA 2x SYBR Green qPCR Mix is considered for high-performance, high-throughput, real-time PCR. The kit covers a Hotstart Taq DNA polymerase engineered through a process of molecular evolution. The result is an exclusive enzyme, specifically designed for qPCR using SYBR Green I dye chemistry.

YTA 2x SYBR Green qPCR Mix is a convenient Mastermix. 2x SYBR Green qPCR Mix is a suitable mastermix of the components (except primers, template, and water) necessary to achieve real-time polymerase chain reaction (PCR) using SYBR Green I dye with enhanced sensitivity and specificity. The SYBR Green dye binds to double-stranded (ds) DNA, thus providing a fluorescent signal that reflects the amount of dsDNA product generated during PCR.

Storage

This reagent can be stored at 4°C for 2 months and protected from light. For longer storage, it should be kept at -20°C and protected from light.

Applications

Gene expression analysis , Low copy gene detection,
Microarray validation ,Gene knockdown validation

Features

- This reagent can be used in glass capillary systems (e.g., LightCycler, Roche Molecular Systems, Inc.).
- Hot Start technology with anti-Taq DNA polymerase antibodies enables high specificity and reproducible amplification.

Composition of the SYBR Green qPCR Mix

100 mM KCl , 5 mM MgCl₂, 400 μM dNTPs, 0.1 U/μl Hot Start Taq DNA Polymerase, 1x SYBR Green and other optimized buffer components.

Protocol

1. Preparation of reaction solution

Add all the solution in a thin walled PCR tube on ice.

For a total 20μl reaction volume

Component of sample	Volume	Final concentration
2x SYBR Green qPCR Mix	10 μl	1X
Forward Primer (10 μM)	0.4 μl	0.2 μM
Reverse Primer(10 μM)	0.4 μl	0.2 μM
Template DNA	variable	(<1ng/μl)
Passive reference Dye(50X) optional-not provided	Variable (following table)	1X
Water, nuclease-free	to 20 μl	–

Passive Reference Dye (50X) Guide Table (Cat No : YT2552/1)	
Instrument	Volume
Rotor-Gene™; DNA Engine Opticon™, Opticon 2, and Chromo 4™ Real-Time Detector; No Rox Mastercycler ep realplex, Smart Cycler , Roche LightCycler 480, Bio-Rad CFX96	No Rox
ABI Step One,ABI Step One Plus, ABI Prism 7000/7300/7700/7900,Eppendorf	0.4 μl
ABI 7500,ABI 7500 Fast, Stratagene Ms3000 / Mx300SP,Rotorgen 3000	0.08 μl

Note:

The primer concentration can be further optimized, if needed. The optimal range for primers is 0.1~1μM. The Template DNA concentration can be further optimized by gradient dilution. The concentration is generally <1ng/μl

2. Setup the plate

- Transfer the appropriate volume of reaction mixture to each well of a PCR tube/plate. Reaction volumes may be scaled down from 20 μl to 10 μl if low volume tubes/plates are used. Cap or seal the reaction tube/plate and centrifuge briefly.

3. Preform qPCR using the following thermal cycling conditions.

Initial Denaturation	95°C	20 sec-3 min*	40 Cycles
Denature	95°C	5 sec **	
Anneal/Extend***	60°C	>20 sec	
Melting curve analysis			

*20 sec at 95°C is sufficient time for enzyme activation, however optimal denaturation of complex

targets may require up to 3 min denaturation.

**Select minimum time (not more than 20 sec) according to instrument user guide.

***For 3 step cycling protocols, anneal at optimal annealing temperature for 20 sec followed by 5 sec extension and data acquisition at 72°C.

4. Analyze the results

Data analysis varies depending on the instrument used. Please refer to your instrument user guide for information.

Important Notes

Template

Genomic DNA, plasmid DNA, or cDNA can be used as template. For optimal quantitative results use up to 20ng of genomic DNA or plasmid DNA per 20µl reaction (for smaller volumes, the amount of template should be decreased equivalently). Using greater amounts of template may reduce the maximum fluorescence signal and linearity of standard curves due to binding of the SYBR Green I dye to the template. For two-step RT-PCR, use either undiluted or diluted cDNA generated from up to 1µg of total RNA. The volume of the cDNA added (from the RT reaction) should not exceed 10% of the final PCR volume (e.g., for a 20µl qPCR reaction, use up to 2.0µl of undiluted cDNA).

Primers

Careful primer design and purification (HPLC-purified primers are recommended) is particularly important in order to minimize loss in sensitivity due to the production of nonspecific amplification products in SYBR Green I-based qPCR. This effect becomes more prominent at low target concentrations. To maximize the sensitivity of the assay, use the lowest concentration of primers that can be used without compromising the efficiency of the PCR reaction (50-400nM of each primer). For optimal results, design primers that amplify PCR products 60-400 bp in length. The primers should exhibit a melting temperature (T_m) of approximately 60°C, to take advantage of two-step cycling. If

performing real-time two-step RT-PCR, we recommend designing primers specifically for amplification of cDNA derived from mRNA. This prevents amplification of contaminating genomic DNA and inaccurate quantification of mRNA.

Melting Curve Analysis

Following real-time qPCR, melting curve analysis should always be performed to identify the presence of primer-dimers and analyze the specificity of the reaction. Program your thermocycler according to the instructions provided.

SYBR Green I

2x SYBR Green qPCR Mix contains an elevated, optimized concentration of the fluorescent dye,

SYBR Green I. High signal intensities are achieved as a result of increased tolerance to high concentrations of SYBR Green I by the engineered, novel SYBR DNA Polymerase. SYBR Green I binds all double-stranded DNA molecules, emitting a fluorescent signal on binding. The excitation and emission maxima of SYBR Green I are at 494 nm and 521 nm, respectively, which are compatible with use on any real-time cyclers.

Magnesium chloride

The MgCl₂ concentration in 2x SYBR Green qPCR Mix is optimized for most primer combinations. You do not need to add additional MgCl₂ to the mix to get efficient and specific PCR.

Guidelines for preventing contamination of qPCR reaction

During qPCR more than 10 million copies of template DNA are generated. Therefore, care must be taken to avoid contamination with other templates and amplicons that may be present in the laboratory environment. General recommendations to lower the risk of contamination are as follows:

- Prepare your DNA sample, set up the qPCR mixture, perform thermal cycling and analyze PCR products in separate areas.
- Set up qPCR mixtures in a laminar flow cabinet equipped with an UV lamp.
- Wear fresh gloves for DNA purification and reaction set up.
- Use reagent containers dedicated for PCR. Use positive displacement pipettes, or use pipette tips with aerosol filters to prepare DNA samples and perform PCR set up.
- Always perform "no template control" (NTC) reactions to check for contamination

Quality Control

The absence of endodeoxyribonucleases, exodeoxyribonucleases and ribonucleases is confirmed by appropriate quality tests. Functionally tested in amplification of a single-copy gene from human genomic DNA.