



Total RNA Extraction mini Kit plus YTzol

Cat No : YT9001

Size: 50 preps

<<for research use only>>

Kit Contents	50 preps
YTzol Lysis buffer	50 ml
Wash buffer 1	30 ml
Wash buffer 2 (concentrate)	20
RNase free water	10 ml
Mini column	50 each
Elution tube	50 each
Collection tube	50 each

Preparation note:

Add 80ml ethanol (96-100) to wash 2

Application:

Real-time-PCR (RT-PCR) , Real-time quantitative,
Northern blotting , Nuclease protection assays ,
RNA amplification for microarray analysis ,
cDNA library preparation after poly (A)+ selection

Materials be supplied by the users : Ethanol (96–100%)

Description : This Special RNA Extraction Kit provides a simple method of isolating total RNA from a wide range of sample types and amounts. In general, samples are lysed and then homogenized in the presence of guanidinium isothiocyanate, a chaotropic salt capable of protecting the RNA from endogenous RNases. After homogenization, ethanol is added to the sample. The sample is then processed through a spin column containing a clear silica-based membrane to which the RNA binds. Any impurities are effectively removed by subsequent washing. The purified total RNA is then eluted in RNase-Free Water and is suitable for use in a variety of downstream applications.

Feature

- Stable yield
- high-quality purified total RNA
- No need for chloroform, phase separation and precipitation steps
- Fast RNA extraction by YTzol

Note

- Use sterile, disposable, and individually wrapped plastic-ware.
- Use only sterile, disposable RNase-free pipet tips and microcentrifuge tubes.
- Wear disposable gloves while handling reagents and RNA samples to prevent RNase contamination from the surface of the skin. Change gloves frequently, particularly as the protocol progresses from crude extracts to more purified material.
- Always use proper microbiological aseptic techniques when working with RNA.
- Perform all extraction steps on ice
- Recommended volume of YTzol

10cm ² adherent cells	1 ml
10 ⁷ suspension cells	1-2 ml
100 ul white cells	2 ml
50-100 mg ordinary tissue	1 ml
50-100 mg special tissue(live,spleen,bone or cartilage)	2ml
15-100 mg plant tissue	1 ml



Protocol

1. Sample process

Tissues from animal or plant(either fresh or frozen at -70°C until use) can be processed by freezing with liquid nitrogen and grinding into a fine powder using a mortar and pestle. Homogenize tissue samples in 500 μL to 1 ml YTzol per 10-50 mg tissue using a tissue homogenizer or rotor-stator.

Adherent Cells : Lyse cells directly in a culture dish by adding 1 ml of YTzol to the dish and passing the cell lysate several times through a pipet tip. The amount of YTzol required is based on the culture dish area (1 ml per 10 cm^2) and not on the number of cells present.

Suspension Cells: Harvest cells and pellet cells by centrifugation. Use 1 ml of the YTzol per $5-10 \times 10^6$ animal, plant, or yeast cells, or per 1×10^7 bacterial cells. Lyse cells by repetitive pipetting up and down. Do not wash cells before addition of YTzol to avoid any mRNA degradation.

Disruption of some yeast and bacterial cells may require a homogenizer.

2. Incubate at $15-30^{\circ}\text{C}$ for 5 min, to lyse the nucleoprotein complex completely.

3. centrifuge or spin ,To remove particulate debris from homogenized tissue and transfer the supernatant into a new nuclease-free tube

4. Add an equal volume of ethanol. Mix well, a visible precipitate may form after adding ethanol.

5. Place a Mini Column to a Collection Tube and transfer the ethanol added sample mixture to the Mini Column. Centrifuge at full speed for 1 min, discard the flow-through and return the Mini Column back to the Collection Tube.

6. Add 500 μL of Wash Buffer 1 to the Mini Column, centrifuge at full speed for 1 min. Discard the flow-through and return the Mini Column back to the Collection Tube.

7. Add 750 μL of Wash Buffer 2 to the Mini Column, centrifuge at full speed for 1 min. Discard the flow-through and return the Mini Column back to the Collection Tube. -- Note: Make sure that ethanol has been added into Wash Buffer 2 when first use.

8. Repeat step 7 for one more washing.

9. Centrifuge the Mini Column at full speed for an additional 3 min to dry the Mini Column. -- Important Step! This step will avoid the residual liquid to inhibit subsequent enzymatic reaction.

10. Place the Mini Column to a Elution Tube (provided, 1.5 ml microcentrifuge tube).

11. Add 20-50 μL of RNase-free ddH₂O to the membrane center of the Mini Column.

Stand the Mini Column for 1 min. -- Important Step! For effective elution, make sure that RNase-free ddH₂O is dispensed on the membrane center and is absorbed completely. -- Important: Do not elute the RNA using RNase-free water less than suggested volume ($< 40 \mu\text{L}$). It will lower the RNA yield.

12. Centrifuge the Mini Column at full speed for 1 min to elute RNA.

13. Store RNA at -70°C .