



### YTA Multiplex PCR MasterMix 2X with UDG

2X Multiplex PCR Mastermix +UDG	100 µl
---------------------------------	--------

**For research use only**

Cat No: YT1557

Size: 1 ml ( 40 rxn)

Store at -20°C or at 4°C for up to 30 days.

YTA Multiplex PCR Master Mix with UDG is used to prevent non-specific PCR amplification and contamination which is containing of all components (except primers and templates) required for multiplex PCR reactions. DUTP /UDG anti-pollution system is introduced in the Mix, and UDG enzyme can rapidly degrade the pollutants containing U at room temperature. Hotstart Taq DNA Polymerase's superior performance in combination with an optimized buffer system increases the specificity of the reaction.

#### Protocol

Note: Before setting up the PCR reactions, prepare a Primer Mix II with 0.5 µM of each primer.

#### Prepare the PCR Reaction Mix

1. Allow all reagents to thaw completely. Mix gently by inverting the tube. Spin briefly. Put all reagents on ice.
2. Using a 96-well optical reaction plate, add the following to one well per sample:

#### 1. Preparation of reaction solution

Add the following reagents to the proper thermal cycler reaction tube or plate on ice:

Component	Volume	Final concentration
2X Multiplex PCR Master Mix with UDG	25 µl	1X
Primer Mix II (0.5 µM each)	5 µl	50 nM each primer[1]
Template DNA	0.1–0.2 µg	2–4 ng/µL
Nuclease-free water	Adjust to 50 µL	n/a

1) The final concentration of each primer in a typical PCR is between 0.05-0.4 µM. In most cases, a final concentration of 0.15 µM gives satisfactory results. Increasing the primer concentration up to 0.4 µM may increase the yield.

2) Seal the reaction plate with Clear Adhesive Film.

#### Amplify DNA for Analysis

Choose an amplification protocol based on your analysis method

#### Amplify for analysis by agarose gel electrophoresis

1. Configure the run method as outlined in your instrument's user manual. Use the following parameters:

Technical whatsapp : +989038623150

Order whatsapp: +989371095037

Office : 021 - 40777399



Stage	Time	Temperature
Hold	3 min	95°C
35 Cycles	30 sec	95°C
	90 sec	58 °C
	90 sec	72 °C
Hold	5 min	72 °C
Hold	□	4

- Mix well and briefly spin the reaction plate.
- Load the reaction plate into the instrument, and start the run. See your instrument’s user manual for detailed instructions on how to load and run the plate.
- Analyze the data according to the instructions in the user manual for your gel electrophoresis instrument.

### Amplify for CfDNA Sequencing

Component	Volume
2X Multiplex PCR Mastermix +UDG	12.5 µl
Primer Mix II	2 µl
cfDNA	X µl
Nuclease-free water	To 25 µl

- Configure the run method as outlined in your instrument’s user manual. Use the following parameters:

Stage	Time	Temperature
Hold	3 min	95°C
25-35 Cycles	30 sec	95°C
	90 sec	58 °C
	90 sec	72 °C
Hold	5 min	60 °C
Hold	□	4

- Mix well and briefly spin the reaction plate.
- Load the reaction plate into the instrument, and start the run. See your instrument’s user manual for detailed instructions on how to load and run the plate.