



cDNA Synthesis Kit

For research use only

Cat No: YT4500

Size: 50 Rxn

Store at -20°C

| Contents: | volume |
|--------------------------------------|-------------------|
| M-MLV (10,000 U) 200u/ μl | 50 μl |
| 5X first-strand Buffer | 200 μl |
| Oligo(dT)18 primer | * |
| Random hexamer primer | ** |
| RNasein (40u/ μl) | 25 μl |
| dNTP 10mM | 50 μl |

Important Note:

Primers are lyophilized.

Please add DEPC water used for Con.50 μM

* Oligo(dT)18 primer : 60.96 μl

** Random Hexamer : 184.48 μl

cDNA synthesis kit offer whole components you need for cDNA synthesise include RNase inhibitor, Primer, dNTP mixture, Buffer and M-MuLV RT enzyme . M-MuLV is a reverse transcriptase enzyme with optimum activity at 37C and no RNase H activity. for greater application flexibility we provide both random hexamer and Oligo dt primers. This kit provides superior reagents that ensure efficient first-strand cDNA synthesis from mRNA or total RNA templates.

Features :

- Have ease of use and consistency
- High sensitivity reverse transcription from low abundance template
- broad linear dynamic range of total input RNA (0.1ng – 5 μg) with a highly efficient MMLV reverse transcriptase
- Complete cDNA synthesis and qPCR on the same day

Specification :

- store at -20°C
- Primers are lyophilized
- Contains both Oligo(dT)18 primer and Random Hexamer
- It can be used with either sequence-specific primers, poly(dT)18primers, or random Hexamer

Applications :

- First-strand cDNA synthesis for RT-PCR and qPCR
- Gene expression profiling
- Reagents can be used for miRNA synthesis

Protocol (50 Test)

First-Strand cDNA Synthesis Using M-MLV RT

1. Add the following reagents into a sterile, nuclease-free nuclease-free tube on ice in the indicated order:

| | | |
|--------------------|--|---|
| Template RNA | Total RNA Poly (A) mRNA Or specific RNA | 0.1ng - 5 μg 1 to 500 ng 1-5 μg |
| Prime | oligo (dT)18 primer(50 μM) or Random hexamer primer(50 μM) | 1.0 μl 1.0 μl |
| DEPC-treated water | | To 13.5 μl |
| Total Volume | | 13.5 μl |

2. Mix gently, centrifuge briefly and incubate at 70°C for 5 min. Chill on ice, spin down and place the vial back on ice.
3. Prepare the following cDNA Synthesis Mix, add the following components in the indicated order:

| | |
|------------------------------|-------------------|
| 5x first-strand buffer | 4 μl |
| dNTPs(10 mM each) | 1 μl |
| RNasin (40U/ μl) | 0.5 μl |
| M-MLV | 1 μl |

4. Mix gently and centrifuge
5. For oligo(dT)18, incubate for 60 min at 42°C . For random hexamer primed synthesis, incubate for 60 min at 37°C .

Technical whatsapp : +989038623150

Order whatsapp: +989371095037

Office : 021 - 40777399



6. . Terminate the reaction by heating at 70°C for 5 min. The reverse transcription reaction product can be used immediately in second strand cDNA synthesis reactions or stored at -20°C for less than a week. For longer storage, -70°C is recommended.

Protocol (100 Test final volume 10µl)

The kit contains high-quality reagents, allowing it to be used for 100 samples. To utilize it, follow these steps:

- 1- Add the following reagents into a sterile, nuclease-free nuclease-free tube on ice in the indicated order:

| | | |
|--------------------|--|--------------------------------------|
| Template RNA | Total RNA Poly (A) mRNA Or specific RNA | 0.1ng - 3µg 1 to 300 ng 1-3 µg |
| primer | oligo (dT)18 primer(50µM) or Random hexamer primer(50µM) | 1µl or 1µl |
| DEPC-treated water | | To 6.75 |
| Total Volume | | 6.75 |

- 2- Mix gently, centrifuge briefly and incubate at 70°C for 5 min. Chill on ice, spin down and place the vial back on ice.
3- Prepare the following cDNA Synthesis Mix, add the following components in the indicated order:

| | |
|------------------------|--------|
| 5x first-strand buffer | 2µl |
| dNTPs(10 mM each) | 0.5µl |
| RNasin (40U/µl) | 0.25µl |
| M-MLV | 0.5µl |

- 4- Mix gently and centrifuge
5- For oligo(dT)18, incubate for 60 min at 42°C. For random hexamer primed synthesis, incubate for 60 min at 37°C .
6-Terminate the reaction by heating at 70°C for 5 min. The reverse transcription reaction product can be used immediately in second strand cDNA synthesis reactions or stored at -20°C for less than a week. For longer storage, -70°C is recommended.
Notice that the final volume of your cDNA will be 10µl.

Protocol (100 Test final volume 20µl)

First-Strand cDNA Synthesis Using M-MLV RT

- 1- Add the following reagents into a sterile, nuclease-free nuclease-free tube on ice in the indicated order:

| | | |
|--------------------|--|--------------------------------------|
| Template RNA | Total RNA Poly (A) mRNA Or specific RNA | 0.1ng - 5µg 1 to 500 ng 1-5 µg |
| Prime | oligo (dT)18 primer(50µM) or Random hexamer primer(50µM) | 1.0 µl 1.0 µl |
| DEPC-treated water | | To 13.5 µl |
| Total Volume | | 13.5 µl |

- 2- Mix gently, centrifuge briefly and incubate at 70°C for 5 min. Chill on ice, spin down and place the vial back on ice.
3- Prepare the following cDNA Synthesis Mix, add the following components in the indicated order:

| | |
|------------------------|-----------|
| 5x first-strand buffer | 2 µl |
| dNTPs(10 mM each) | 0.5 µl |
| RNasin (40U/µl) | 0.25 µl |
| M-MLV | 0.5 µl |
| Nuclease free water | Up to 6.5 |

- 4- Mix gently and centrifuge
5- For oligo(dT)18, incubate for 60 min at 42°C. For random hexamer primed synthesis, incubate for 60 min at 37°C .
6- Terminate the reaction by heating at 70°C for 5 min. The reverse transcription reaction product can be used immediately in second strand cDNA synthesis reactions or stored at -20°C for less than a week. For longer storage, -70°C is recommended.