



## Tissue Genomic DNA Extraction mini Kit

(Cat: FATGK001)  
(For research use only)

### Special protocol for extraction of DNA from milk sample

#### Important Notes:

1. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
  2. Check FATG1 Buffer before use, Warm fATG1 Buffer at 60°C for 10 minutes if any precipitate formed.
  3. Add indicated volume of ethanol (96-100%) to Wash Buffer before use.
  4. Prepare a heating block or a water bath to 60 °C.
  5. All centrifuge steps are done at full speed (~18,000 x g) in a microcentrifuge.
  6. Preheat Elution Buffer or ddH<sub>2</sub>O to 65°C for elution step.
  7. If the sample is containing few DNA and getting more yield is needed, please preheat the Elution Buffer.
- We recommend 65°C as the heating temperature.

#### Additional requirement:

Lysozyme  
Benchtop microcentrifuge  
Vortex or bead beater equipment  
1.5 ml microcentrifuge tubes  
Heating block or a water bath to 60°C and 70°C  
96~100 % ethanol

#### Please Read Important Notes of FATGK User Manual Before Starting the Following steps.

Hint: • Set dry or water baths: 37 °C for step 2 and 60 °C for step 4.

- Preheat Elution Buffer to 65 °C for elution step

**1A.** Transfer **up to 1 ml of milk sample** to a microcentrifuge tube (not provided) and centrifuge at full speed (~18,000 x g) for 3 minutes. Discard the supernatant including the creamy floats on the top layer after centrifugation and use a paper towel or a cotton swap to remove any white remains on the tube wall.

**1B.** Transfer **up to 20 mg of powder milk sample** to a microcentrifuge tube (not provided). Add 1 ml of ddH<sub>2</sub>O and mix well by vortexing. Centrifuge at full speed (~18,000 x g) for 3 minutes. Discard the supernatant including the creamy floats on the top layer after centrifugation and use a paper towel or a cotton swap to remove any white remains on the tube wall.

2. Add **450 µl of Elution Buffer and Lysozyme (not provided) to final 0.5 mg/ml**. Resuspend the pellet and mix the mixture by vortexing. Incubate at **37°C for 30 min**.

3. Centrifuge the sample at full speed for 3 min and discard the supernatant completely by pipetting.

4. Add **450 µl of FATG1 and 45 µl of Proteinase K solution (10 mg/ml)**. Resuspend the pellet and mix the mixture by vortexing. Incubate at **60 °C for 30 min**.

5. Add **450 µl of FATG2** and mix the sample mixture by vortexing. Incubate at **70 °C for 15 min**.

6. Add **450 µl ethanol (96~100%)** to the sample mixture. Mix thoroughly by pulse-vortexing for 5 sec.

7. Place a FATG Column to a Collection Tube.

8. Transfer **up to 750 µl** of the sample mixture to the FATG Column and centrifuge at full speed for 1 min. Discard the flow-through and place the FATG Column back to the Collection Tube.

9. Repeat Step 8 for the rest of the sample mixture.

10. Add 400 µl W1 Buffer to the FATG Mini Column. Centrifuge at full speed for 1 min then discard flow-through.

---Make sure that ethanol has been added into W1 Buffer when first open.

11. Add 650 µl Wash Buffer to the FATG Mini Column. Centrifuge at full speed for 1 min then discard flow-through.



## Tissue Genomic DNA Extraction mini Kit

- Make sure that ethanol has been added into Wash Buffer when first open.
12. Centrifuge at full speed for an additional 3 min to dry the column.
- **Important Step! This step will remove the residual liquid.**
13. Add 40~100  $\mu$ l of preheated Elution Buffer or ddH<sub>2</sub>O (pH 7.5-9.0) to the membrane of the FATG Mini Column. Stand the FATG Mini Column for 3 min.
- **Important Step!** For effective elution, make sure that the elution solution is dispensed onto the membrane center and is absorbed completely.
- If less sample to be used, reduce the elution volume to 40  $\mu$ l to increase DNA concentration and do not elute the DNA using less than suggested volume (40  $\mu$ l). It will lower the final yield.
14. Centrifuge at full speed for 2 min to elute DNA.