



## DEPC ( Diethyl pyrocarbonate )

## For research use only

Cat No : YT9079

Size : 5 ml

Store at 2 – 8 °c

### How to make DEPC-treated water :

To make DEPC (Diethyl Pyrocarbonate)-treated water, which is used to create RNase-free water for RNA work, follow these steps:

### Materials Needed:

- DEPC (Diethyl Pyrocarbonate)
- Distilled water (dH<sub>2</sub>O)
- Glass container
- Magnetic stirrer or orbital shaker
- Autoclave

### Procedure:

1. **Preparation:**
  - Add **1 mL of DEPC** to **1 L of distilled water** in a glass container for a final concentration of 0.1% DEPC
  - Mix well to ensure proper dispersion of DEPC in the water.
2. **Incubation:**
  - Allow the mixture to react at **room temperature for at least 1 hour** or at **37°C for 12 hours**.
  - Stir the solution using a magnetic stirrer or orbital shaker to ensure thorough mixing.
3. **Autoclaving:**
  - Autoclave the DEPC-treated water at **15 psi** for 15-45 minutes ( **20 minutes is recommended**) to inactivate any remaining DEPC.
  - Allow the solution to cool to room temperature after autoclaving.
4. **Storage:**
  - Store the DEPC-treated water at room temperature in a clean, aseptic container.
  - Label the container clearly to avoid confusion with untreated water.

### Safety Note:

- DEPC is a strong RNase inhibitor but can be harmful if inhaled or ingested. Always handle it with care, use appropriate personal protective equipment (PPE), and work in a well-ventilated area.

### Applications of DEPC-Treated Water

1. RNA Isolation
2. cDNA Synthesis
3. PCR and qPCR
4. RNA Sequencing (RNA-Seq)
5. In Vitro Transcription
6. Northern Blotting
7. RNA Stability Studies