

dNTP Mix 10 mM

Cat No : YT5553 Size : 0.4 ml

Cat No : YT5554 Size: 0.4 ml X5

YTA dNTP Mix (10mM) is created by enzymatic synthesis techniques and superior purification protocols ensure the highest quality for deoxynucleotides. This product is ultrapure (>99%) and the quality is checked by a set of PCR, RT-PCR, and Klenow reactions. YTA dNTP Mix is suitable for various PCR amplifications including long PCR.

Storage temperature:

Stable for 2 years in the case of storage at -20, or 7 days at 25°C before activity begins to reduce.

For research use only